

ROMA



Power Pack

Model No. **P001**

Special Features

- Twin drive wheels for greater traction
- Dismantles easily for easy transport and storage
- Easy to operate hand control lever
- Variable speed setting with reverse
- Security key switch
- Removeable battery pack for remote charging
- Lifts out of use to enable normal wheelchair function
- Suitable for users up to 136 kgs / 21 stone



Specifications



Model No.*	P001	
Battery Capacity	12v / 20Ah	
Charger (off board)	1.5A	
Drive Wheels (puncture proof)	2 x 200 mm / 8"	
Maximum Drive Speed	3.4 mph / 5.5 kmph	
Maximum Range	Up to 10 miles / 16 kms	
Fits Wheelchair Width	360-510 mm / 14-20"	
Overall Weight	14 kgs / 31 lbs	
Motor Weight	6.4 kgs / 14 lbs	
Maximum User Weight	136 kg / 300 lbs**	

* Lower limb amputees should refer to user manual before use. ** Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.